

# CT Myocardial ECV Analysis

Extracellular Fluid Volume display allows evaluation of myocardial fibrosis

Non-contrast and delayed acquisition images are used to display subtracted myocardium and calculate ECV. Color maps are shown with polar maps for endo- and epicardial regions.

## Key features:

- Select multiple data for analysis using non-contrast and contrast exams
- Workflow and automated extraction of myocardial walls gives quick results
- Editing tools for adjusting automatically extracted axis and contours
- Apply hematocrit value to calculate ECV

