

# Multiparametric MR Analysis

## High-precision time intensity curves with the use of deformable registration

Using contrasted dynamic phase T1 MR images, parametric maps can be displayed, and curves are shown in general time intensity or as defined by the BI-RADS category. Color maps showing wash-in and washout lets you easily see locations of lesions. Images of the same position from multiple phases can be viewed at the same time. With deformable registration, you can see the time intensity curves with higher precision.

### Key features:

- Load from 3 to as many 40 phases from a single series
- Automatically assigns pre-contrast, early and late phase images
- Wash-in and washout display options
- Time-intensity curves generated with deformable registration
- Graph axes can display time or phase, and intensity or percent
- MIP MPR window generates a reference image
- Exportable analysis and key images for BI-RADS compliant reporting

